



*stoma-belt. The new belt for greater security and freedom.*

*Safe. Even when travelling.*

*stoma-belt is supplied in a high-quality travel case. Unobtrusive and practical, there is also room for all the other supplies you need for stoma care.*

*The stoma-belt is available in the following versions:*

*Velcro fastener / black: □ 129 incl. 16% VAT.  
Sizes: XS, S, M, L, XL, XXL, XXXL*

*Hook fastener / black: □ 129 incl. 16% VAT.  
Sizes: S, M, L, XL*

*Dispatch is cash on delivery and comes to you entirely free of charge.*



**stoma-belt**



*stoma-belt is available from salomé-rökona, rökona Textilwerk GmbH,  
Schaffhausenstr.101, 72072 Tübingen, Germany  
Telephone: +49 (0) 7071 / 153-404, Fax: +49 (0) 7071 / 3 52 41.  
Also online under [www.stoma-belt.com](http://www.stoma-belt.com)*



**stoma-belt**

## *Confidence.*

*stoma-belt is entirely new, developed by patients for patients under the leadership of Achim Krüger, himself a stoma patient. Which is why stoma-belt guarantees problem-free use in **any situation**.*

***Clinically tested and successfully tried** by stoma patients, stoma-belt lets you feel **totally secure** and gives you unlimited **physical freedom** as no other belt can.*

*For a new  
lease of life.*

## Security.

**stoma-belt gives you all-round reliable protection and a sense of absolute security.**

Experience renewed **security and comfort** with stoma-belt. The **elasticated material** ensures that the pouch **fits snugly and securely** in the stoma-belt and can therefore be worn **almost invisibly** under clothing.

It can be emptied and the **pouch safely changed** in just two easy steps.

The breathable, **moisture-absorbent material** provides a high level of wearing comfort and **stops odour**. Both day and night.

## Closeness.

**stoma-belt is exceptionally discreet and hygienic.**

Because it is extremely **hard-wearing**, stoma-belt adjusts to every body shape and every situation, guaranteeing **maximum compression** on the body without bulkiness or rustling.

Even extreme movements, such as when running or swimming, or also intimate physical contact present no barrier. Discover a **new lease of life** with stoma-belt.

## Wake up to new life.

**stoma-belt. Confidence. Security. Closeness.**

- Clinically tested
- Highly elasticated material
- Flexible and adjustable to every abdominal size
- No bulkiness or rustling under clothing
- Flange presses perfectly onto the abdominal wall
- Individual adjustment and positioning for one and two-piece pouch systems
- Pouch fits snugly and securely in the high-elasticity integrated pocket
- No build-up of odour
- Waterproof
- Easy-care
- Compatible with all commercially available pouch systems


**stoma-belt is available in black with Velcro fastener in sizes XS, S, M, L, XL, XXL, XXXL; with hook fastening in sizes S, M, L, XL.**

stoma-belt is available from salomé-rökona.


For further information, please contact our hotline: +49 (0) 7071 / 153-404. Achim Krüger will be waiting to give you detailed advice.

## Simple. Practical. Secure.

stoma-belt is supplied with a hole already provided for the retaining flange or other flange-pouch links. If required this hole may easily be enlarged by cutting along the coloured seams with scissors. Please note that the outer black seam must not be altered.



**1)** With the pouch already fixed to the body, place the stoma-belt over the pouch with the fastener to the front. Close the Velcro or hook fastener and turn the stoma-belt until the pocket is positioned over the pouch.



**2)** Guide the pouch into the belt pocket through the custom-cut hole.



**3)** Insert one hand through the lower (zipped) slit of the pocket and insert the other hand through the upper slit and pull the pouch into shape until it lies flat in the pocket.



**4)** Adjust the pouch to the length of the pocket by rolling up the narrow end, and then close the zip.



**5)** stoma-belt will now press the skin barrier firmly to your body and the pouch will lie securely fixed in the belt pocket.